Distributed by:

(352) 377-1140 □ (866) 284-9690 □ www.exac.com

Processed by:

(386) 418-8888 □ (800) 624-7238 □ www.rtx.com

REGENAFIL® and REGENAFORM® ALLOGRAFT PASTES

DESCRIPTION

All paste products contain donated human tissue for transplantation. Flowable allograft paste contains human demineralized bone matrix (DBM) in an inert porcine gelatin carrier. Moldable allograft paste contains human DBM and human cortical cancellous chips (CCC) in an inert porcine gelatin carrier.

These products were evaluated in a human clinical study and were shown to induce bone formation. Each lot of product is tested using the athymic nude rat assay to verify osteoinductive potential.

Caution: United States of America Federal law restricts this device to sale by or on the order of a licensed physician, dentist, or practitioner.

Table 1: Composition of Products

PRODUCT	DBM	CCC	Gelatin	Water
REGENAFIL® Allograft Paste	X		X	X
REGENAFIL® RT Allograft Paste	X		X	
REGENAFORM® Allograft Paste	X	X	X	X
REGENAFORM® RT Allograft Paste	X	X	X	

DONOR SCREENING AND TESTING

This donated human tissue was recovered using aseptic technique from a donor determined to be eligible for tissue donation. After consent for donation is obtained, tissue donors are screened for risk factors associated with infectious diseases and medical conditions that would rule-out donation. This screening includes, but is not limited to: Family/next-of-kin interview, medical/hospital record evaluation, donor physical assessment and autopsy report (if performed). The final determination of donor eligibility is made by a licensed physician utilizing all available, relevant information.

Donor medical charts are reviewed by an RTI Biologics (RTI) Medical Director for donor eligibility. Following processing, production records are reviewed by qualified RTI personnel and the tissue determined to be suitable for transplantation. The donor's blood is tested by a CLIA certified laboratory utilizing test kits approved by the FDA for donor screening and cadaveric blood testing where applicable. The following test criteria were met for this donor:

Table 2: Required infectious disease tests and acceptable result

BLOOD TEST	Acceptable Result
HIV-1 / HIV-2 Antibody	Non-Reactive/Negative
Hepatitis C Virus Antibody	Non-Reactive/Negative
Hepatitis B Surface Antigen	Non-Reactive/Negative
Hepatitis B Core Antibody (Total)	Non-Reactive/Negative
RPR for Syphilis	Non-Reactive/Negative
Human T-Cell Lymphotropic Virus I/II Antibody	Non-Reactive/Negative
HIV-1/HCV NAT-TMA (tested after March 2005)	Non-Reactive/Negative
HIV-1 DNA by PCR (tested prior to March 2005)	Non-Reactive/Negative

PROCESSING

Tissues are processed into their final forms in controlled environments. Microbial testing is performed where appropriate and results meet a documented acceptance criterion that is based on AATB and FDA requirements.

STERILE R PRODUCTS

Products terminally sterilized with gamma irradiation are treated with a validated* dose.

*All references to "Validated" sterilization processes indicate that the validation meets or exceeds requirements for product sterilization, i.e., a SAL of 10⁻⁶ per AAMI and ISO standards.

ADVERSE REACTIONS

Surgeons should discuss these potential adverse events with their patients.

The following complications may occur after tissue transplantation:

- Loss of function and/or integrity of transplanted tissue due to resorption, fragmentation, and/or disintegration including, but not limited to, associated loss of continuity or displacement.
- Transmission of known infectious agents including, but not limited to, HIV, Hepatitis B, Hepatitis C, syphilis and bacteria.

S00722

- Immune rejection of transplanted tissues; or allergic reactions to residual chemicals.

This product may contain trace amount of any of the following chemicals: ascorbic acid, hydrochloric acid, hydrogen peroxide, hydroxymethyl aminomethane, isopropyl alcohol, povidone-iodine, sodium phosphate buffer.

The same medical/surgical conditions that may complicate any surgical procedure may occur during or following the transplantation of this product.

Exactech and RTI, their affiliates, distributors and agents disclaim all implied warranties concerning this tissue including any warranty of merchantability or fitness for particular purpose. Please promptly report complaints and possible adverse events as instructed in the Customer Service section of this insert.

INDICATIONS FOR USE

These products are intended to be packed into bony voids or gaps to fill and/or augment dental intraosseous, oral and cranio-/maxillofacial defects. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone, including Periodontal defects; Alveolar ridge augmentation; Extraction sockets (ridge preservation, implant preparation/ placement); Maxillary sinus floor elevation; craniofacial augmentation; Root resection, apicoectomy and cystectomy; Tumor resection.

One or more of the product formulations, depending upon specific anatomical location and physician and/or dentist preference, can be placed in the dental intraosseous defect site.

CONTRAINDICATIONS

These products are contraindicated where stabilization of the defect is not possible.

WARNINGS

A small number of patients may experience localized immunological reactions resulting from the use of this product that generally consist of transient localized edema and swelling, and rash. Although there is no evidence that the product will be unsafe or ineffective in such patients, the safety and effectiveness of the product in these patients has not been established. These products are contraindicated, and should not be used, in patients who have known allergies to porcine collagen or who are being treated for desensitization for meat products, as this product contains porcine collagen.

These products do not possess sufficient mechanical strength to support the reduction of a defect site prior to soft and hard tissue ingrowths.

RECEIPT INSPECTION

Upon arrival, unpack and inspect the graft. Verify that you have received the correct graft type and that the expiration date has not elapsed. Should you receive the incorrect order or find that the shipping container and product packaging integrity is compromised, do not use the product and immediately notify Exactech at (866) 284-9690 or (352) 377-1140.

STORAGE

Storage conditions are also stated on the package label. It is the responsibility of the end-user to maintain products intended for transplantation at appropriate storage conditions. Storage conditions should be documented.

Table 3: Storage of Products

Product Type (FZ / DR)	Temperature	Expiration
Frozen (FZ) on label	-40 to -85°C (-40 to -121°F)	See the expiration date on the label
Frozen (FZ) on label	-20 to -39°C (-4 to -39°F)	Either 6 months or the expiration date on the label, whichever is shorter
Room Temperature (DR) on label	1 to 37°C (33 to 98.6°F)	See the expiration date on the label

GENERAL INFORMATION

- Once opened, this allograft paste must be used for the current procedure or discarded.
- This allograft paste may not be re-sterilized.
- Do not use expired or damaged product.
- Tissue & package contents for single patient & single surgery use only.
- Adhere to instructions provided on surgical techniques and use of instrumentation.
- Additional products of varying types and sizes should be accessible in case of an unexpected need.
- This product and all packaging materials used are latex free.
- The box or clear plastic outside packaging acts as a dust cover and is used to protect the product during shipping and storage. Outside dust cover of frozen product may be vented or perforated.
- The surgeon should be familiar with the graft, surgical procedure and any instruments prior to use of this graft or subsequent products.

PRODUCT INFORMATION AND PREPARATION:

- After thawing or mixing paste with warm water, DO NOT DIRECTLY expose paste to water warmer than 38°C (100°F). This may cause the paste to dissolve. Frozen paste must be thawed inside the sealed inner foil pouch.
- Remove the double packaged allograft product, package insert and Tissue Utilization Record card.
- Compare the labeling identification number on the pouch with the number on the outer pouch. If the numbers do not match, do not use the tissue.
- Use sterile technique to peel open the outer package and pass the inner package containing the product and any accessories to the sterile field for further handling.

DIRECTIONS FOR USE AND MIXING:

FROZEN FLOWABLE AND MOLDABLE ALLOGRAFT PASTE

These products must be warmed prior to surgery. Carefully follow the preparation steps described below.

- If the paste is packaged in a dispenser, open the inner package in the sterile field, and lay out the contents. If the paste is NOT packaged in a dispenser, DO NOT OPEN the inner package; the paste must be warmed while still sealed.
- Prepare either a standard warmer with sterile drape or sterile thermos. **Thermos:** Add a solution of sterile water or sterile saline that is between 43 and 49°C (110 to 120°F). **Warmer:** Warm the sterile water/saline between 43 and 49°C (110 to 120°F). **NOTE:** DO NOT USE WATER WARMER THAN 49°C (120°F). (DO NOT heat any frozen paste for more than six hours).
- Place the capped dispenser or sealed package into thermos or warmer.
- Wait approximately 5 minutes for the paste to warm; then remove.
- Uncap or open package and use as needed.
- If the paste does not extrude or mold, it needs further warming. Recap dispenser or place paste into watertight sterile container and warm for 5 minutes.

ROOM TEMPERATURE FLOWABLE ALLOGRAFT PASTE

REGENAFIL® RT ALLOGRAFT PASTE can be mixed with sterile saline, sterile water or the patient's own blood. The following table indicates the recommended volume of fluid required to mix.

Table 4: Approximate Fluid Volumes for Mixing

Room Temperature Flowable Allograft Paste	
Product & Fluid Volumes in (cc)	0.5cc
Whole blood (Patient's blood)	Red Line (0.5)
Room temperature (19 to 25°C; 66 to 77°F)	Black line (0.3)
Warm (43 to 49°C; 110 to 120°F)	Black line (0.3)

Room Temperature flowable allograft paste must be mixed with (warm or room temperature) fluids or patient's blood prior to implantation.

- Open the package containing the fluid syringe and lay out the contents on the sterile field. Repeat for the powder syringe.
- Pick up powder syringe; do not remove cap. Draw back the plunger of the powder syringe about 3/4 of complete extension. Then shake the powder syringe vigorous manner for about 10 seconds. This ensures proper mix of product that may have settled during storage.
- To reduce the trapped air in the powder, gently push in the powder syringe stem to compress the powder; careful not to compact powder in excess. Compacted powder can reduce ability of proper mixing.
- Fill the fluid syringe with room temperature (19 to 25°C; 66 to 77°F) sterile saline or sterile water to black line. If using the patient's blood, fill to the red line. Tap the syringe to eliminate air bubbles.
- NOTE:** (For more robust handling use warm (43 to 49°C; 110 to 120°F) sterile fluids; fill to the black line. Fluids may be warmed in a warmer or thermos.)
- Remove cap from powder syringe and hold in vertical position. Holding the uncapped syringe with the tip down will result in loss of product. Connect the fluid syringe with the powder syringe keeping the fluid syringe vertical and above the powder syringe.
- Push the fluid from the fluid syringe slowly into the powder syringe. Apply counter pressure to the plunger of the powder syringe to keep the powder compressed while transferring the fluid. Fluid should penetrate the powder. Allow fluid to reach base of powder; it is okay for some fluid to remain in fluid syringe before mixing.
- Hold for 10 seconds, and then turn syringes to a horizontal position for mixing.
- Push mixture back and forth between syringes, empty each side as much as possible during each stroke. Continue mixing for 10 cycles (back and forth equals one cycle). Paste mixture should all transfer from one syringe to the other once fully mixed.
- Push all product into the fluid syringe.
- Connect preferred tip, and dispense as needed.

ROOM TEMPERATURE MOLDABLE PASTES

REGENAFORM® RT ALLOGRAFT PASTE can be mixed with sterile saline, sterile water or the patient's own blood. The following table indicates the approximate volume of fluid required to mix with RT moldable allograft paste:

Table 5: Approximate Fluid Volumes for Mixing

Room Temperature Moldable Allograft Paste		
Product & Fluid Volumes in (cc)	1cc	2cc
Whole blood & Room Temperature (19 to 25°C; 66 to 77°F)	Red line (0.6)	Red line (1.2)
Warm (43 to 49°C; 110 to 120°F)	Blue line (0.4)	Blue line (0.9)

Room temperature moldable allograft paste must be mixed with (warm or room temperature) fluids or patient's blood prior to implantation.

- Open the inner package and place the contents onto the sterile field.
- Shake package containing powder, then open.
- Fill the fluid syringe with room temperature (19 to 25°C; 66 to 77°F) sterile saline or sterile water to red line. If using the patient's blood, fill to the red line. Tap the syringe to eliminate air bubbles.

NOTE: (For more robust handling use warm (43 to 49°C; 110 to 120°F) sterile fluids; fill to the blue line. Fluids may be warmed in a warmer or thermos.)

- Add the fluid into the powder slowly.
- Use spatula to mix powder and fluid completely.
- Paste is now ready for molding.
- Mold graft into the desired shape (e.g., disc or log). At room temperature, the graft should become firm.

NOTE: Moistening the surgical glove will prevent paste adhesion to the glove.

TISSUE UTILIZATION RECORD (TUR CARD)

The TUR card is designed for tracking implanted grafts. All information provided is kept confidential and used for product tracking only. This card must be completed by the operating room staff or clinician after the surgery. For further instructions please see those listed on the TUR card. Return completed copy to RTI to address or fax provided on TUR card.

FREQUENTLY ASKED QUESTIONS

Q. If we don't use all of the paste at once, can we use it later in the same case?

A. Yes. Replace the syringe cap on the syringe, tighten securely, and place it in sterile warm water at 43 to 49°C (110 to 120°F). Alternatively the product may be warmed by placing it in a sealed watertight sterile container and then into sterile warm water at 43 to 49°C (110 to 120°F).

Q. Can over-warming the paste damage it?

A. Yes. As with all DBM pastes, exposing the paste to temperatures in excess of 54°C (129°F) for over 30 minutes will denature some of the proteins and destroy the product.

Q. What if the paste stops coming out of the syringe or stops molding easily mid-use?

A. The paste may have cooled. Recap the syringe (if not packaged in a syringe, place in a watertight sterile container) and warm for approximately 5 minutes (or until paste flows or shapes easily) before reusing.

Q. Are frozen and room temperature allograft paste products tested for osteoinductivity?

A. Yes. A sample from every lot is implanted in vivo in an athymic nude rat to test its osteoinductive potential. Only lots that pass this quality control test are released for distribution. Findings from an animal model are not necessarily predictive of human clinical results.

Q. What does the FZ or DR mean on the label?

A. FZ stands for frozen and appears on all frozen products. DR stands for dried or freeze dried and refers to room temperature products.

CUSTOMER SERVICE

For further information or to report a complaint or adverse event, please contact:

Exactech, Inc.
Gainesville, Florida, 32653, USA.
TEL: (352) 377-1140, (866) 284-9690 or FAX: (800)-329-3922
www.exac.com

Table 6: Definitions of label symbols:

 See instructions for use	 Sterile by Gamma Irradiation	 Expiration Date (month/day/year)
 Catalog number	 Lot number (Donor number)	 Serial Number (Tissue number)
 Storage temperature limits	 Do not reuse	 Manufacturer

REGENAFIL® and REGENAFORM® are registered trademarks of RTI Biologics, Inc.

 RTI Biologics, Inc.
11621 Research Circle | Alachua, FL 32615-2650 U.S.A
www.rtx.com

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E Exactech	RTI BIOLOGICS™
(352) 377-1140 • (866) 284-9690 • www.exac.com	+1 386.418.8888 • +1 800.624.7238 www.rti.x.com

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REGENAFIL® e REGENAFORM® sono marchi depositati di RTI Biologics, Inc.