


Cartilage Processor

Reveille™



Surgeon focused. Patient driven.™

Exactech®
B I O L O G I C S

BENEFITS/FEATURES:


- Precision crafted blade and sieve design permits high speed cutting
- Designed to create convective movement of fluid during processing
- Easy to assemble; with detachable tissue reservoir
- Two minute processing time yields tissue particles with increased surface area suitable for grafting procedures
- Fluid environment protects tissue viability without the need for cell culture
- Easy to load and retrieve tissue
- High cell viability
- Produces tissue to fill defects up to 4cm²
 - Host tissue eliminates the potential of graft rejection
 - Off-the-shelf convenience
 - Increase in surface area allows greater chondrocyte exposure




Crafted Blade for High Speed Cutting



Precision in Size Reduction and Separation
Particle Size Approximately 800 to 1200 Microns¹



35x: Underneath Surface Layer



1500x: Underneath Surface Layer

Increase in Chondrocyte Exposure

REFERENCES

1. Data on File in Exactech Taiwan.

Exactech is proud to have offices and distributors around the globe.

For more information about Exactech products available in your country, please visit www.exac.com.

Reveille[™] Cartilage Processor is manufactured by Exactech Taiwan, Ltd, and distributed by Exactech, Inc.