

DESIGN TEAM

- **Constantine Demetrapoulos, MD**
*Hospital for Special Surgery
New York, N.Y.*
- **James K. DeOrio, MD**
*Duke University
Durham, NC*
- **Mark E. Easley, MD**
*Duke University
Durham, NC*
- **James A. Nunley, MD**
*Duke University
Durham, NC*
- **Jason Nowak, DPM, FACFAS**
*Shasta Orthopaedics
Redding, CA*
- **Robert Santrock, MD**
*Orthopaedic Foot & Ankle Center
Worthington, Ohio*

SCOPE & MATERIALS

- Tibia and Talar Cutting Guides
- Fluoroscopic Insert
- Reciprocating Saw Blade
- Case Report
- Anatomic Model

TARGETING

- Fellowship-trained foot and ankle MDs (special focus on those using PROPHECY™ Navigation)
- Current Vantage Ankle users
- Foot and ankle MDs
- DPM's certified in rear-foot surgeries and currently performing total ankle replacements

SELLING FEATURES

- 3D-printed nylon cutting guides designed to uniquely fit each patient
- Delivered within approximately three weeks*
- Pre-plan tibia and talar resections to align and cut each in one step
- Used in combination with the Vantage Total Ankle System
- Designed to provide OR efficiency, reliability and aid soft tissue preservation
- Designed for a reproducible fit

LIST PRICING

- Tibia and Talar Cutting Guides **\$2,000**

ORDERING INFORMATION

To upload your uncompressed DICOM file format scans, **visit www.vantagepsi.com**. Please ensure all scans are uploaded a minimum of 21 days before surgery.

Refer to the Vantage Ankle PSI protocol (00-0000616) for complete CT scan/X-ray details.

Ensure you place an order for standard Vantage Ankle Instrumentation and implants (350A, 350B, 351A, 351B and 351F)

CATALOG #	DESCRIPTION
351-50-00	Vantage Ankle PSI Instrument Kit

COMPETITION

Wright PROPHECY Pre-Operative Navigation System

- Used with INFINITY™, INBONE™, INVISION™

Paragon 28 MAVEN™ Patient-Specific Instrumentation

- Used with APEX 3D™ Total Ankle Replacement System

**Orders will be fulfilled within three weeks as long as there are no extraordinary circumstances.*

