

# A Solution for Every Glenoid Challenge.

**equinox<sup>®</sup>**

Platform Shoulder System



## You're Prepared. So Are We. A Solution for Every Challenge.

Our newest implants are engineered to solve clinical challenges and improve patient outcomes.



### POSTERIOR AUGMENT CAGE GLENOID (8 DEGREE)

- Initial fixation through press-fit bone cage
- Preserves bone by building up the eroded side rather than reaming to correct the deficiency
- Allow for fixation of the true axis of the glenoid



### SUPERIOR/POSTERIOR AUGMENT BASEPLATE

- Corrects glenoid deformities with both superior and posterior wear
- Preserves bone by building up the eroded sides rather than reaming

A wide range of augmented glenoid options is available, offering bone preserving solutions for challenging glenoids.



**10 Degree Superior Augment Baseplate**



**Extended Cage Baseplate**



**8 Degree Posterior Augment Baseplate**



**8 or 16 Degree Posterior Augment Glenoid**



**Cage Glenoid**



**+4mm Expanded Glenospheres**

Exactech is proud to have offices and distributors around the globe. For more information about Exactech products available in your country, please visit [www.exac.com](http://www.exac.com).