

equinox[®]
PRESERVE STEM



exactech


equinox[®] | PRESERVE STEM

The Equinox[®] Preserve stem is created from the collaborative efforts of our global surgeon thought leaders and world-class engineers and designed to preserve distal humeral bone while also achieving optimal stability and fixation. This platform stem is compatible with the Equinox Shoulder System that uniquely focuses on anatomical replication, reverse shoulder arthroplasty and options for difficult glenoids and revisions.

REVISION FRIENDLY

- Flexibility to choose between anatomic and reverse.
- Streamlined technique with fewer steps
- Shortened stem length and polished distal tip in case of stem removal

Plasma coating designed for press-fit fixation



Grit blast coating to allow for initial fixation



Tapered stem to facilitate insertion

Full-length fins are designed to provide rotational stability.¹



PRESERVE AND SMALL REVERSE

Courtesy of Jeffrey Hartzell, MD

REFERENCES: 1. Data on File at Exactech, Inc.